

## Catalog Number

40031-PG6

**Guard Column starter set, UDA (wcx)**, complete with 5 each 4.0 mm ID x 10 mm long cartridges packed with 4  $\mu$ m, 100 A phase and 1 each guard column holder. Cogent TYPE-C Brand. 1 PK.



### Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, delrin case to protect and keep it from drying out. Use with [guard column holder](#) sold separately.

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

*Use this guard column cartridge with a 3.0, 4.0, and 4.6 mm ID analytical columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$241.85**

Printed On 2024/05/18