

## Catalog Number

40200-15P-2

**HPLC Column, Silica-C 2.0**, 2.2um, 2.1mm ID x 150mm Length, 120A.  
Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Cogent Silica-C 2.0™ “near UHPLC”, HPLC Column is made with 2.2um, 120A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

*This column can be used for LCMS when increased resolution is desired.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,503.14**

Printed On 2024/05/31