


Catalog Number

45099-110

Polymer Tubing Drill. 1.0mm Polymeric Pin, Accessory for Flanging Tool. ARE-Applied Research brand. 1 EA.



Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool and is not included in the Kit or with the Tool.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.68

Printed On 2024/05/20