

## Catalog Number

45099-115

**Polymer Tubing Drill. 1.5mm Polymeric Pin, Accessory for Flanging Tool.** ARE-Applied Research brand. 1 EA.



### Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool and is not included with the tool or with the kit.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$45.68**

Printed On 2024/05/31