

Catalog Number

46000-HP

Tool Kit for HPLC. High Pressure. Includes Tools, Tubing, Fittings, Ferrules, Filters and box. ARE-Applied Research brand. 1 EA.



Additional Info

This convenient kit puts all the parts tools you need to keep an HPLC up and running. A must have for everyone that is responsible for maintenance of the instruments.

Includes:

- 10 Ferrules, PEEK, Double. 1/16"
- 10 Ferrules, SS, Single, 1/16"
- 6 Tubing Connector, One Piece, Finger Tight
- 1 Inline Filter, PEEK complete. 2 um PEEK frit.
- 4 Tubing Connectors, Finger Tight Double Cone Ferrules 10-32
- 4 Tubing Connectors, Hex Head, Long for Double Ferrule 10-32
- 5 Tubing Connectors, SS, Long Head, Rheodyne Type, 10-32
- 5 Tubing Connectors, SS, Extra Long, Rheodyne Type, 10-32
- 2 Solvent Filters, Last Drop, SS Element, 2 um

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$879.84

Printed On 2024/05/20