

## Catalog Number

46010-BT

**Beta Tool.** For 1/16" Stainless Steel HPLC Tubing, Red Grips. MicroSolv Brand. 1 EA.



### Additional Info

This is a very handy tool for every laboratory that works with 1/16" Stainless Steel Tubing. These Tubing Pliers / Tubing Straightener has soft grips and a Center Hole that is slightly larger than 1/16" that is designed to remove kinks or bends in tubing or for making seamless bends when you squeeze the Pliers. This will help to insure a good internal flow path in your HPLC tubing and avoid problems due to bent or kinked tubing.

The Beta Tool™ can also be used as precision pliers.

Click [HERE](#) for a snub nose version for smaller spaces.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$176.38**

Printed On 2024/05/20