


## Catalog Number

47305-25-30

**Tubing, PEEK, Pre-Cut, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD, Ultra-T grade, blue intermittent stripes, 30 cm long (+/- 1 mm). ARE-Applied Research brand. 5 PK.**



### Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut PEEK Tubes are the perfect solution.

This Ultra-T™ grade (*Premium*), PEEK tubing is extruded from virgin PEEK material without adding any additives. The broken yellow stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color. The internal and outer dimensional tolerances of this tubing is +/- 0.025mm.

- Color Code: Blue, Stripes for 0.010" (0.25mm) ID.
- Wall Thickness: 0.675 mm
- Pressure rating: 5600 psi.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$72.23**

Printed On 2024/05/19