


## Catalog Number

47305-30-C

**Tubing, PEEK, 0.030 inch (0.75 mm) ID, 1/16th inch (1.6 mm) OD,**  
HPLC grade, green stripe, 25 meter roll, ARE-Applied Research brand. 1  
EA.



### Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The green stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.425 mm
- Color code: green stripe
- Pressure rating: 3500 psi (241.3 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$481.50**

Printed On 2024/05/20