


Catalog Number

47305-40-1C

Tubing, PEEK, 0.040 inch (1.0 mm) ID, 1/16th inch (1.6 mm) OD,
HPLC grade, grey stripe, Your custom continuous length tubing, priced
per meter. ARE-Applied Research brand.



Additional Info

You determine the length of this HPLC grade tubing that is extruded from virgin PEEK material without any manufacturing additives. The grey stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.30 mm
- Color code: grey stripe
- Pressure rating: 2400 psi (165.5 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/30