


## Catalog Number

47305-40-C

**Tubing, PEEK, 0.040 inch (1.0 mm) ID, 1/16th inch (1.6 mm) OD,**  
HPLC grade, grey stripe, 25 meter roll, ARE-Applied Research brand. 1  
EA.



### Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The grey stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.30 mm
- Color code: grey stripe
- Pressure rating: 2400 psi (165.5 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$585.00**

Printed On 2024/05/19