


## Catalog Number

47305-ST-02-1C

**Tubing, PEEK, 0.0025 inch (0.064 mm) ID, 1/16th inch (1.6 mm) OD,** Ultra-T grade, natural color, Your custom continuous length tubing, priced per meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this Ultra-T™ grade (*premium*) tubing that is extruded from virgin PEEK material without any manufacturing additives. A natural PEEK material has no color leachables.

- Choose between a 1 & 300 meter continuous length.
- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.768 mm
- Color code: Natural (no stripe)
- Pressure rating: 6300 psi (434.4 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

*This tubing is available through our Customer Service department only, see link below. Cutting charge applied under 3 meter lengths.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/31