


Catalog Number

47305-ST-05-C

Tubing, PEEK, 0.005 inch (0.13 mm) ID, 1/16th inch (1.6 mm) OD,
Ultra-T grade, red intermittent stripes, 25 meter roll, ARE-Applied
Research brand. 1 EA.



Additional Info

This Ultra-T™ grade (*premium*) tubing is extruded from virgin PEEK material without any manufacturing additives. The broken red stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.735 mm
- Color code: Red stripes
- Pressure rating: 6100 psi (420.6 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$770.40

Printed On 2024/05/20