


## Catalog Number

47310-05-1C

**Tubing, Capillary, PEEK, 75 um ID, 360 um OD, natural color, Your custom continuous length tubing, priced per meter. ARE-Applied Research brand.**



### Additional Info

You determine the length of this small OD, PEEK, capillary tubing. It is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.157 um
- Color code: natural.
- Pressure rating: 5200 psi (358.5 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/20