

Catalog Number

47310-05-1M

Tubing, Capillary, PEEK, 75 um ID, 360 um OD, natural color, 1 meter roll. ARE-Applied Research brand. 1 EA.



Additional Info

This small OD, PEEK, capillary tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Color code: natural
- Pressure rating: 5200 psi (358.5 bar)

Other lengths available:

- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$32.90

Printed On 2024/05/20