


## Catalog Number

47310-05-5M

**Tubing, Capillary, PEEK, 75 um ID, 360 um OD, natural color, 5 meter roll.** ARE-Applied Research brand. 1 EA.



### Additional Info

This small OD, PEEK, capillary tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.157 um
- Color code: natural
- Pressure rating: 5200 psi (358.5 bar)

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$144.45**

Printed On 2024/05/20