

Catalog Number

47310-13-3M

Tubing, PEEK, 0.005 inch (0.13 mm) ID x 1/32nd inch (0.79 mm) OD,
general grade, natural, 3 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives. No color is added leaving the tubing all-natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.33 mm
- Color code: Natural
- Pressure rating: 5500 psi (380 bar)

Other lengths available:

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$104.00

Printed On 2024/05/20