

Catalog Number

47310-18-1C

Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/32 inch (0.79 mm) OD,
general grade, natural, Your continuous length tubing. ARE-Applied
Research brand.



Additional Info

You determine the length of this general grade tubing that is extruded from virgin PEEK material without any manufacturing additives. No color is added leaving the tubing all-natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.305 mm
- Color code: Natural
- Pressure rating: 5100 psi (350 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$34.67

Printed On 2024/05/20