

## Catalog Number

47310-75-1M

**Tubing, PEEK, 0.030 inch (0.75 mm) ID, 1/8th inch (3.2 mm) OD,**  
HPLC grade, natural color, 1 meter roll, ARE-Applied Research brand. 1  
EA.



### Additional Info

This wide OD, PEEK, tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 1.225 mm
- Color code: natural
- Pressure rating: 5000 psi (344.7 bar)
- *Caution Note: This tubing has very thick walls and can unwind very quickly. Handle with care to prevent injury. Safety glasses and gloves should be worn when unpacking.*

Other lengths available:

- [3 meter roll](#)
- [5 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$46.39**

Printed On 2024/05/20