

Catalog Number

47452-14

Adapter, for use with low pressure. Female luer lock to barbed connection for 3/32 inch tubing. Polypropylene. ARE-Applied Research brand. Bulk Price.



Additional Info

This adapter has a female Luer Lock port on one side and a hose barb connector for 3/32-inch (2.4 mm) ID tubing on the other side. *Bulk packaging.*

This adapter is excellent when using a syringe and you want to connect to 3/32-inch tubing.

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for economical pack of 10.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/20