

Catalog Number

49027-04W

End Plugs, for HPLC Columns. Nylon 10-32 Thread. White. ARE-
Applied Research brand. 10/EA.



Additional Info

Use these nylon, endstopper fittings to close and seal any guard, analytical or prep HPLC column with 10-32 threads and prevent it from drying out or from solvent leaking from the column. *Compatible with any brand of column.*

Use when storing columns or when columns are not attached to instruments will increase the column lifetime.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$8.66

Printed On 2024/05/19