

HSQ Kit.

49101-07

Part No.

HSQ Kit, Complete. 75mm HPLC column, test solutions, accessories & software for performance qualification. Chemical Solutions brand. 1 EA.



Additional Info

This complete *high-speed* qualification (HSQ™) kit for self-compliance of HPLC performance qualification (PQ) provides everything you need to complete a full qualification of your HPLC instruments and perform holistic and modular AIQ tests - in less than half a workday. You can do preventative maintenance with immediate instrument uptime without waiting for vendors to verify your instrument. Fulfill your requirements to “maintain instruments in a qualified state.” See USP <1058>. You can use this kit to qualify your HPLC Instruments for precision and accuracy.

The HSQ Kit includes:

- 6 Linearity Solutions - 30 ml (L2 and L3 60 ml)
- 1 Resolution Test Mixture - 30 ml
- 1 Wavelength Calibration Solution (Holmium Oxide) - 4 ml
- 1 Gradient Visualization Solution - 30 ml

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/20