

## Catalog Number

49115-01

**UV-vis Calibration Test Kit, Complete.** Includes NIST traceable solutions, software. Chemical Solutions brand. 1 EA.



### Additional Info

This UV-vis Calibration PQ Kit™ is a series of NIST traceable solutions with easy data recording and reporting of the actual data with acceptance criteria. This will allow you to qualify a UV-vis spectrophotometer easily and fully in regulated laboratories. In-house performance qualification (PQ) of these instruments is easy and saves money and time.

UV-vis Kit includes:

- Absorbance Solutions - 7 x 30 ml each
- Reference Solution Diluent - 60 ml
- Wavelength Calibration Solution - 30 ml
- Stray Light Solution - 30 ml
- Software: Excel Template for Data Recording and Analysis
- Certificate of Analysis.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,250.06**

Printed On 2024/05/31