

Catalog Number

49115-25

UV-vis Calibration Solutions. Replacement set of solutions for use with UV-vis calibration kit. Chemical Solutions brand. 1 EA.



Additional Info

This set of NIST traceable solutions are replacements for use with the Chemical Solutions UV-vis Calibration Kit™ for performance qualification (PQ) of UV or UV-vis spectrophotometers and detectors.

Set includes:

- Absorbance Solutions - 7 x 30 ml each
- Reference Solution Diluent - 60 ml
- Wavelength Calibration Solution - 30 ml
- Stray Light Solution - 30 ml
- Certificate of Analysis.

Click [HERE](#) for more information about the solutions and the kit or Click [HERE](#) for specifications and concentrations.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,241.49

Printed On 2024/05/20