

## Catalog Number

49313-508SE

**Tubing, PEEK, Coiled, 0.020 (0.50 mm) ID, 1/16th inch (1.6 mm) OD,**  
6930 mm length, HPLC grade, orange stripe, with tubing connectors.  
ARE-Applied Research brand. 1 EA.



### Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 250 mm wide when resting and 1300 mm when fully extended, with a total tubing length of 6930 mm, end to end.

*Type 10 FlexiTube, color coded: Orange striped for 0.020 inch (0.50 mm) ID*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$123.26**

Printed On 2024/05/20