

## Catalog Number

49330-F2

**Frit, PEEK, 2  $\mu$ m**, encased in a PEEK ring, 0.20 inch OD .102 inch ID, 0.062 inch thick, for use with 2.1mm HPLC columns. EPP-Porous Products Brand, 1. EA.



### Additional Info

These all PEEK frits are made for use for the ARE-Applied Research brand, all PEEK, 2.1 mm ID HPLC column hardware.

Dimensions are: *0.20 inch (1.57 mm) thick and 0.102 inch (2.59 mm) diameter. These frits have a 2  $\mu$ m porosity and the PEEK ring retains biocompatibility.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$40.45**

Printed On 2024/05/19