

Catalog Number

49330-GC1

Guard Column Kit, PEEK, 3.0 mm ID, 8 mm active length. ARE-Applied Research brand. 1 EA.



Additional Info

This Guard Column Hardware Kit includes one 3.0mm ID Column Tube, two frits and the Guard Column Holder. All components, including the Frits are extruded from ® PEEK. No chemical treatment is done to the tubes and the interior mirror like finish is provided through our very careful extrusion techniques.

Carbon is added to black components.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$126.90

Printed On 2024/05/19