


## Slurry Reservoir R10

49331-SR1  
Part No.

**Slurry Reservoir Kit, All-PEEK, for packing 4.6 mm ID columns that are less than or equal to 100 mm in length.** Reservoir is 100 mm long with a 7.5 mm ID, complete with all parts and accessories needed to



### Additional Info

This kit is for packing All-PEEK columns 4.6 mm ID that are up to 100 mm in length.

The Reservoir Body included in this kit is 100 mm long with 7.5 mm ID and has an internal volume of approximately 4.4 ml. This reservoir is equipped with 5/8-18 external threads on both sides. The adapter has an internal 1/2-20 screw thread on one side and internal 5/8-18 screw threads on the other side. The adapter connects the reservoir to the column being packed.

Kit includes everything you need to pack a bioinert All-PEEK column:

- One Slurry Reservoir - the "body" which is black PEEK.
- One [End Fitting](#) to connect tubing from a pump to the reservoir.
- Two [PEEK sealing ring / frit spacer](#) for the bottom and top connection.
- One [Packing Adapter](#) to connect the reservoir to the column.
- One [PEEK sealing ring](#) for the top of the column connection.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$857.08**

Printed On 2024/05/20