

## Catalog Number

49430-100

**Tubing, Stainless Steel, Pre-Cut**, 0.030 inch (0.762 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, White Band, 100cm Long. ARE-Applied Research brand. 1 EA.



## Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut Stainless Steel Tubes are the perfect solution.

Welded and annealed, the ends are cut using Electronic Discharge (*EDC*) Technique to achieve square and burr free zero dead volume connections. They are tested according to EN 10204/3.1B. The White Band is a color code for .030" ID.

- *Internal and Outer Dimension Tolerances:* +/- 0.05mm.
- Wall Thickness: 0.414 mm

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$202.98**

Printed On 2024/05/19