


## Catalog Number

58610-PT47

**Filters, Membrane, 47mm, 1.5um, PTFE.** Round, MicroSolv Brand.  
10/PK.



### Additional Info

These PTFE, MicroSolv Brand, Filter Discs can be used with the Joshua Direct Filter Degasser or any device when membrane filtration is desired.

Hydrophobic or water resistant membranes with minimal extractables and absolute retention, PTFE is chemically resistant and will not degrade with Normal Phase Solvents. They can be autoclaved at 121°C for 30 minutes or EO sterilization without damage to the integrity or performance of the filter.

- Diameter: 47mm
- Max Operating Temperature: 90°C
- Pore Size: 1.5 um

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$136.60**

Printed On 2024/05/19