

Catalog Number

58710-S-C

Syringe, Chrom. Disposable, 10ml, Luer Slip, Non Sterile. MicroSolv Brand. 100/PK.



Additional Info

These Luer Slip Tip, *Chrom Syringes* are chosen for many uses in chromatography as well as used with syringe filters or for any lab situation needing an inert, non-reactive, plastic syringe.

Latex free, they do not contain rubber, silicone oil, styrene or DEHP that is found in other plastic syringes. These all plastic syringes are manufactured from laboratory grade polypropylene and polyethylene. More chemically resistant than syringes with rubber-tipped plungers, they are also equipped with a positive "click", safety stop to help prevent accidental spills.

These economical, disposable syringes are not individually wrapped and are bulk packed in one bag.

Click [HERE](#) for smaller package quantity or Click [HERE](#) for case quantity.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$84.29

Printed On 2024/05/20