

## Catalog Number

58802-01

**Bottle, Mobile Phase Reservoir, Safety, Coated Glass, Clear, 1L, Screw Top, GL45 Threads. MicroSolv Brand. 1 EA.**



### Additional Info

This Clear, Plastic Coated Borosilicate Glass Bottle has a 1L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, one liter Capacity and White Enamel Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#), other Degassers or disposable [Filter Degassers](#). A Solid Screw Cap without a liner and a drip ring is included. This Bottle is often used with the [Mobile Phase Delivery System](#).

*The outside only of the Bottle is coated in a transparent, thick epoxy plastic and if this bottle is dropped and broken, the plastic will prevent a spill incident and prevent workers from getting cut with glass or splashed with solvent.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$0.00**

Printed On 2024/05/20