


## Catalog Number

58802-30

**Pouring Ring**, Replacement, Blue, for HPLC Mobile Phase Reservoir Bottles. GL45mm Threads. ARE-Applied Research brand. 1 EA.



### Additional Info

Using this Polypropylene Pouring or "Drip Ring" can make decanting liquids from any bottle with a GL45 Thread easier and prevent that last drop from getting away from you.

Friction fit makes this ring stay in the bottle during pouring and capping.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1.93**

Printed On 2024/05/20