

## Catalog Number

58807-100-A-CS

**Bottles, HPLC Mobile Phase Storage**, Glass, Amber, 100ml, Screw Top, GL45 Threads. MicroSolv Brand. 10/CS.



### Additional Info

These Amber, Borosilicate Glass Bottles have a 100ml "Total Fill Capacity" and are perfect as HPLC Mobile Phase Storage Bottles. The small openings and Enamel Graduations and Marking Patches makes them a convenient way to store concentrated Mobile Phase or its Ingredients.

Blue, Solid Screw Caps without liners and polypropylene Drip Rings are included. This Bottle is *Not "Safety Coated"*.

Click [HERE](#) for smaller package quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$938.18**

Printed On 2024/05/20