

## Catalog Number

58807-2L-CS

**Bottles, Mobile Phase Reservoir**, Glass, Clear, 2L, Screw Top, GL45 Threads. MicroSolv Brand. 5/CS.



### Additional Info

Some of our Most Popular Bottles, these Clear, Borosilicate Glass Bottles have a 2L "Total Fill Capacity" and are perfect as HPLC Mobile Phase Reservoirs. The small openings, one liter Capacity and White Enamel Graduations on the outside of the Bottles help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#), other Degassers or disposable [Filter Degassers](#). A Solid Screw Cap without a liner and a polypropylene Drip Ring are included. This Bottle is often used with the [Mobile Phase Delivery System](#) and is *Not "Safety Coated"*.

Click [HERE](#) for smaller package quantity.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$299.74**

Printed On 2024/05/20