

Catalog Number

58807-45

Pouring Ring, Replacement, Natural, for HPLC Mobile Phase Reservoir Bottles. GL45 Threads. MicroSolv Brand. 1 EA.



Additional Info

Using this Natural color Polypropylene Pouring or "Drip Ring" can make decanting liquids from any bottle with a GL45 Thread easier and prevent that last drop from getting away from you.

This ring will fit deeper into the mouth of the Bottle, making is more secure and less likely to fall out.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2.26

Printed On 2024/05/20