

Catalog Number

58819-05R

Bottle, Mobile Phase Reservoir, Safety, Coated Glass, Clear, 2L,
Screw Top, GL45 Threads, Open Red Cap. MicroSolv Brand. 1 EA.



Additional Info

This Clear, Plastic Coated, Borosilicate Glass Bottle has a 2L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, two-liter Capacity and White Enamel Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#), other Degassers, disposable [Filter Degassers](#). A Red, Open Top Screw Cap and a Red polypropylene Pouring Ring are included.

This Bottle is often used with the [Mobile Phase Delivery System](#). The outside only of the Bottle is coated in a transparent, thick epoxy plastic and if this bottle is dropped and broken, the plastic will prevent a spill incident and prevent workers from getting cut with glass or splashed with solvent.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$255.53

Printed On 2024/05/20