

Catalog Number

59050-25P

HPLC Column, MicroBee Amino, 10um, 4.6mm ID x 250mm Length,
125A. Cogent Brand. 1 EA.



Additional Info

This Cogent MicroBee Amino™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is not end capped and manufactured with traditional bonding of a 3 carbon chain ligand with a terminal Amino group, to irregular shaped silica. This Column is an ideal Alternate in USP Methods that use the L8 Classification and Microbondapak™ or similar Amino columns.

This Column is packed in and shipped in Normal Phase Solvents.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$764.59

Printed On 2024/05/20