

Catalog Number

59051-25P

HPLC Column, MicroBee Cyano, 10um, 4.6mm ID x 250mm Length, 125A. Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent MicroBee Cyano™ Stationary Phase is part of the Similar™ product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A", 10um, 125A Stationary Phase is not end capped and is manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond. Bonded to irregular shaped silica, this makes it an ideal Alternate in USP Methods that use the L10 Classification and Microbondapak™ or similar Cyano Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$764.59

Printed On 2024/05/19