

Catalog Number

60022-05P

HPLC Column, RP PFP, 5um, 4.6mm ID x 50mm Length, 60A. Cogent Brand. 1 EA.



Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines; the highest purity silica available is used making this column ideal as an Alternate in USP Methods that require the L43 Classification.

With a 5µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Often chosen for samples with halogenated or polar functional groups.

Packed and shipped in Reversed Phase Solvents.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$738.25

Printed On 2024/05/20