

Catalog Number

62508-05P-2

HPLC Column, Hyper Similar C8 Octyl, 5um, 2.1mm ID x 50mm
Length, 150A. Cogent Brand. 1 EA.



Additional Info

This Cogent Hyper Similar C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L7 Classification and Thermo Fisher Hypersil™ Columns or Similar C8 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids but expect less retention than with a C18 Column.

A very popular column size, this short and narrow-bore column provides short run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes. This column size is excellent, and a favorite for LCMS.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$785.53

Printed On 2024/05/19