

Catalog Number

62518-25P

HPLC Column, Hyper Similars C18, 5um, 4.6mm ID x 250mm Length, 150A. Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent Hyper Similars C18™ Stationary Phase is part of the Similars product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L1 Classification and Thermo Fisher Hypersil™ Columns or Similar C18 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

This is one of our Most Popular Hyper Similars Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/05/19