

## Catalog Number

63198-C2

**HPLC Column, C18 ODS 3, 10um, 4.6mm ID x 250mm Length, 100A**  
with WVS Cartridge hardware. Partisil Brand. 1 EA.



### Additional Info

This Partisil C18™ Octadecyl Silane (ODS3), Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is fully end capped and modified with traditional bonding of a C18 ligand. This Column is written in many USP Methods; Partisil ODS3 used in older methods with the L1 Classification.

*Often used for hydrophobic or non polar compounds in Reversed Phase HPLC, this column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires end fitting kit [66900-WEF](#).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,399.75**

Printed On 2024/05/20