

Catalog Number

64418-30P-3.9

HPLC Column, Nova Similars C18, 4um, 3.9mm ID x 300mm Length.
Cogent Brand. 1 EA



Additional Info

This Cogent Nova Similars C18™ Stationary Phase is part of the Similars product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very Similar surface area to Nova-Pak®. This makes this column an ideal Alternate in USP Methods that use the L1 Classification and Waters Nova-Pak™ Columns or Similar C18 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,041.49

Printed On 2024/05/20