

Catalog Number

64424-25P

HPLC Column, Nova Similars Phenyl, 4um, 4.6mm ID x 250mm
Length. Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent Nova Similars Phenyl™ Stationary Phase is part of the Similars product line that has been a premium value used by HPLC laboratories for many years.

This "Type-B" Stationary Phase is fully end capped along with bonding of a 4 carbon ligand with a terminal Phenyl group to spherical shaped silica that has very Similar surface area to the original Nova-Pak®. This makes this column an ideal Alternate in USP Methods that use the L11 Classification and Waters Nova-Pak™ Columns or Similar Phenyl Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/05/31