

Catalog Number

65308-25P-4

HPLC Column, Zx Similar C8 Octyl, 3um, 4.0mm ID x 250mm Length.
Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent Zx Similar C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to spherical shaped silica with very Similar surface area to Agilent's Zorbax® columns. This makes it an ideal Alternate in USP Methods that use the L7 Classification and Agilent's Zorbax™ brand C8 HPLC Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids with less Retention than a C18 Column.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,036.93

Printed On 2024/05/20