

## Catalog Number

65518-25P

**HPLC Column, Zx Similar C18, 5um, 4.6mm ID x 250mm Length.**

Cogent Brand. 1 EA.

**MOST POPULAR.**



### Additional Info

This Cogent Zx Similar C18™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very Similar surface area to Zorbax® columns. This makes it an ideal Alternate in USP Methods that use the L1 Classification and Agilent's Zorbax™ brand C18 HPLC Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

*This is one of our Most Popular Zx Similar Columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$921.70**

Printed On 2024/05/19