

## Catalog Number

66086-C1

**HPLC Column, SCX**, 5um, 4.6mm ID x 125mm Length, 100A with WVS  
Cartridge hardware. Partisphere Brand. 1 EA.



### Additional Info

This Partisphere SCX™ Strong Cation Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped and is modified with a ligand that contains sulfonic acid. This Column is written in many USP Methods that call for a Partisil SCX Column and the L22 Classification.

*Often used for purification or analysis of biological compounds in Ion Exchange HPLC, this column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires the end fitting kit [66900-WEF](#) which is sold separately.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,060.03**

Printed On 2024/05/19