

## Catalog Number

67086-15P-4

**HPLC Column, IEX SCX**, 10um, 4.0mm ID x 150mm Length, 100A.  
Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent Ion Exchange Column is a Strong Cation Exchanger (SCX has a negatively charged surface) Stationary Phase that has been a premium value and used by HPLC laboratories for many years. It is manufactured with Sulfonic Acid functional groups that are always charged.

We use only the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L9 Classification or for "Bases" analysis or purification.

With a 10µm Particle, this column works well with higher flow rates and lower back pressures for challenging IEX applications that include Cations (positive charge compounds).

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$785.53**

Printed On 2024/05/20