

## Catalog Number

68004-05P-4-3M

**HPLC Column, RP C4 Butyl, Wide Pore**, 10um, 4.0mm ID x 50mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available. This makes this column ideal as an Alternate in USP Methods that use the L26 Classification, for large molecules.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for proteins, peptides, enzymes and other large, hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for bioanalysis.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$576.05**

Printed On 2024/05/20