

## Catalog Number

68012-P21-PG1

**Prep Guard Column Cartridge, RP C12, 10um, 100A. 21.2mm ID x 10mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Guard Column Cartridge is packed with 10µm, Cogent RP C12™ Stationary Phase. It is double end capped (*fully*), with bonding of a rigid, Adamantine shaped carbon silane (diamond like) to the highest purity silica available. It is designed to protect your RP C12 Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

*Use this Guard Column Cartridge with a 21.2mm ID Prep Column Holder. Click [HERE](#) for more information.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/19